

# Genius Total Solution 8080

Номер части: XX-XXXX



## Описание

Genius 1053 и Genius 3070 построены на основе проверенной технологии Peaks, обеспечивая комплексное газовое решение для тройного квадрупольного Shimadzu LCMS-8080.



## Generator Details

- **Типы газа:** Азот
- **Встроенный компрессор:** Нет

## Ключевая особенность

**Гибкость** - доступен с внутренним воздушным компрессором или без него. **Чистота** - нулевой воздух с содержанием углеводородов менее 0,1 ppm (в виде метана). **Долговечность** - мало движущихся частей, требующих минимального обслуживания. **Решение из одного источника** - один генератор может удовлетворить требования нескольких приложений. **Простая установка** - генератор выполнен в виде plug & amp; игровая система **Экономичный** - более экономичный, чем любой другой способ подачи газа. **Удобно** - газ по запросу, без опасности для здоровья, не нужно беспокоиться о том, что бензин закончится.

## Технические характеристики

### Genius Total Solution 8080

Тип газа	Азот
Время запуска:	60 mins
Текущий:	0
Уровень шума:	0

# [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

