

# ZA070A Zero Air

Номер части: UK: 30-1610 (230v) | EU:  
30-1620 (230v) | US: 30-1630 (110v)




## Описание

Генерируя нулевой воздух для детекторного газа в приложениях ГХ, генератор ZA070A удаляет все углеводороды из источника сжатого воздуха до того, как он попадет в ваше оборудование, обеспечивая самые чистые результаты от вашего детектора.



## Generator Details

- **Типы газа:** Воздух
- **Встроенный компрессор:** Нет



## Ключевая особенность

- Обеспечивает постоянный поток воздуха нулевого качества с содержанием углеводородов менее 0,1 частей на миллион.
- При необходимости работает круглосуточно и без выходных.
- Один генератор может питать несколько приборов.
- Меньшее количество движущихся частей означает минимальное техническое обслуживание и обеспечивает долгий срок службы генератора.
- Требуется минимальная настройка.
- 12 месяцев полной гарантии.

## Технические характеристики

	ZA070A
Тип газа	Воздух
Макс. расход газа:	7L/min
Содержание углеводородов:	0.1ppm
Максимальная чистота:	0.1ppm THC
Время запуска:	40 Minutes
Напряжение:	230 ± 10% VAC
Частота:	50 HZ
Текущий:	4.9 Amps
Максимальная рабочая температура:	30°C / 86°F
Частицы:	0.01µm
Размер (ВхШхГ), мм	880 x 430 x 410 mm
Размер (ВхШхГ), дюймы	34,6 x 16,9 x 16,1 inches
Вес генератора	64kg / 140,8lbs

# [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

